# Ever Simpler:


Streamlining Proteomics with S-Trap TM Turbo and BCA-No-More TM

Sandra Wilson<sup>1</sup>; Darryl J.C. Pappin<sup>1,2</sup>; Alexandre Zougman<sup>2</sup>; John Wilson<sup>1</sup>

<sup>1</sup>ProtiFi LLC, Fairport, NY, United States; <sup>2</sup>University of Leeds, Leeds, United Kingdom




ECENT ADVANCES IN PROTEOMICS enable high-throughput sample identification and quantification, necessitating improved sample preparation workflows. To enhance throughput and robustness, simplifying workflows is crucial. The widely used S-Trap™ system has standardized proteomics sample preparation. Building on this, we introduce a streamlined workflow with S-Trap™ Turbo and BCA-No-More™ technologies to eliminate steps of protein quantification and concentration (Fig. 1C and E) to further simplify and speed up sample preparation.



# Workflow Simplification-

-TRAP™ PLATES AND COLUMNS EFFICIENTLY CAPTURE proteins while removing vexatious contaminants like buffers, salts, reducing and capping reagents, detergents, polymers and small molecules. These impurities can interfere with protein assays, digestion, and MS analysis.

-TRAP™ TURBO IMPROVES UPON THE ORIGINAL design with polymeric materials featuring over 100 times greater derivatization density (Fig. 2A to B). This enhances binding efficiency, reduces the volume of captured matrix required, and yields small elution volumes which eliminats the need for additional concentration steps. Elutions are compatible with direct LC-MS injection.




CA-NO-MORETM INTEGRATES PROTEIN quantification directly into the same plate used for downstream processing, eliminating both the need for a separate assay and sample loss. This time-saving method accommodates a wide range of lysis buffers that would otherwise disrupt sample preparation. Proteins loaded onto S-Traps™ are initially captured at the top of the column (Fig. 3), becoming more concentrated as the uppermost affinity sites fill. This surface-concentration of


Fig. 3 proteins with intrinsically fluorescent tryptophan, tyrosine, and phenylalanine residues (C) allows fluorescent protein quantification through top excitation (A) and emission (B) detection: more protein results in more

## Methods -

-TRAPTM 96-WELL PLATES and new snap-cap S-Trap™ Turbo micro columns followed standard protocols for lysis, reduction, alkylation, denaturation, binding, washing, and tryptic digestion. Various sample hydrophobicities were tested including: serum




HE OPTIMAL Z-POSITION (2100 μm) for protein concentration quantification via fluorescence on the S-Trap™ 96-well plate was determined using a Tecan Spark plate reader. Protein fluorescence was measured in wet and dry states, with excitation between 269 and 280 nm and emission between 325 and 475 nm in top-read mode.




# Results

0

FTER DIGESTION, over 80% of the total peptide amount was recovered in the first elution fraction, equivalent to roughly 60% of the initially applied protein mass (Fig. 6). The columns exhibited high binding efficiency, capturing approximately 90% of applied protein.



ECOVERY: S-TRAPTM TURBO COLUMN: Recovery remained steady regardless of increases in elution volume from 5 μL to 10 μL or 25 μL (Fig. 7) or protein amounts from 1 µg to 100 µg (Fig. 6 and 8). In terms of digestion completeness, peptide yield, identification numbers and detected proteins, S-Trap<sup>TM</sup> Turbos yielded statistically identical results compared to standard S-Traps<sup>TM</sup>.



ECOND AND THIRD ELUTIONS yielded an additional ~10% of the total applied protein (Fig. 8). Consistent elutions as low as 5 μL were achieved (Fig. 7). Recovery was consistent across sample types, volumes did not necessitate speed-vaccing and S-Trap™ Turbo elutions could

be directly loaded onto an autosampler with or without acidification.


#### LUORESCENT CONDITIONS: Optimal

excitation and emission were found to be 277 nm and 350 nm, respectively, except for very low protein amounts where 410 nm was found to yield the best signal to noise. BCA-No-More™ quantification exhibited a linear response up to approximately 100 µg per well. Beyond this, the response decreased as


fluorescent moieties bound deeper within the UV-opaque trapping matrix, forming a curve with a "hook" (Fig. 9).

**UANTIFICATION:** Protein quantification showed strong

microdissection, and an amount not accessible via colorimetric assays. Absorbative losses of such limited amounts were mitigated by handling proteins in 5% SDS during



PROTEIN [µg]



## Conclusion +

## A STREAMLINED, HIGH-THROUGHPUT WORKFLOW

Bringing precision omics to life!

https://protifi.com | info@protifi.com

HIS COMBINED APPROACH PROMISES A SIGNIFICANT reduction in the time and money from sample receipt to analysis readiness, crucial aspects for advancing high-throughput proteomic research and applications.


### SIMPLIFIES WORKFLOW



LIMINATES REDUNDANT STEPS: This method removes two steps of the sample preparation process by integrating the parts into one seamless workflow, and one consumable.

EDUCES PROCESSING TIME: The BCA assay steps – sampling, transferring, incubating, and reading – are eliminated with S-Trap™ Turbo. We avoid time-consuming and -variable sample concentration, greatly shortening preparation and analysis time.

#### ENHANCES EFFICIENCY & ACCURACY



EDUCES EXPERIMENTAL ERROR: This integrated method eliminates errors from sample aliquoting for protein quantification. Simply adjusting the inject volume to match on-column protein loads further decreases potential errors and processing time. This streamlined process minimizes variability and enhances throughput and accuracy.

NHANCES REPRODUCIBILITY: Simplifying procedures, minimizing steps, and reducing manual handling lead to more consistent and reliable


#### SAVE SAMPLE & MONEY



O SAMPLE SACRIFICE, LESS ABSORBATIVE LOSS: The non-destructive nature of the BCA-No-More™ approach preserves precious sample material to enable sensitive quantification at or below the single µg level from any buffer. Additionally, concentration is eliminated and absorbative losses are minimized. This capability is particularly important for low-level analyses and allows distinction between instrument or inject errors and missing initial sample.

## OST-SAVINGS:

Simplifying sample preparation, eliminating and automating steps reduce time, cost, hands-on work and consumable consumption. The risk of costly experimental errors (and repeats) is minimized, further cutting expenses (Fig. 12).



correlation and similar accuracy between traditional BCA assays ( $CV \sim 5\%$ ,  $R^2 \approx 0.99$ , Fig. 11A) and on-plate BCA-No-More<sup>TM</sup>  $(CV \sim 3\% - 5\%, R^2 \approx 0.99, Fig. 11B)$ . This method enabled detection of protein concentrations as low as 1 µg (Fig. 10), an amount for commonly encountered in

samples like in laser capture

 $R^2 = 0.98773$ 

μg PROTEIN

6.0%

2000

intense emission.