

Recommendations for Different Proteases

S-Trap's™ work with your choice of protease.

Protease	Protease:Protein (w/w)	Digestion Time	Digestion Temperature	Digestion Buffer	Digestion Buffer for Isobaric Labeling
Lys-C ¹	~1:30	4 hours	37 °C	Tris-HCl Buffer, pH 8.5	50 mM TEAB
Asp-N ²	1:20	3 hours	37 °C	0.5 mM ZnCl ₂ in 50 mM NH ₄ HCO ₃	0.5 mM ZnCl ₂ in 50 mM TEAB
Glu-C ³	1:10	16 hours	37 °C	0.1 M Phosphate Buffer, pH 7.7*	
Chymotrypsin ⁴	1:30	16-18 hours	37 °C	50 mM TEAB, 1 mM CaCl ₂	
PNGase F (with trypsin) ⁵	1:20	16-18 hours	37 °C	100 mM TEAB, 1 mM CaCl ₂	

*Sodium phosphate must be removed prior to MS-analysis, especially when performing phospho-enrichment

Standard ProtiFi conditions are robust and different samples and experimental variations, especially targeted assays, may benefit from optimization.

Keep in mind:

1. If maximal enzyme/substrate performance is desired, for example to obtain maximum liberation of a target MRM peptide, hold the sample type and amount constant, and test different enzyme:substrate ratios (e.g. 1:10, 1:50, and 1:100 at different times (e.g. 1- 3 hours).
2. A typical protease:protein substrate ratio is 1:20 (w/w). A range of 1:10-1:100 can be implemented.
3. A typical incubation period at 37 °C is 4- 16 hours (overnight). A temperature of 47 °C for a shorter duration (1- 2 hours) can also be used, enzyme dependent.
4. Specifics for your process must be tested accordingly.
5. Proteases:
 - a. Lys-C (Reference: Glatter, T., et al. Large-Scale Quantitative Assessment of Different In-Solution Protein Digestion Protocols Reveals Superior Cleavage Efficiency of Tandem Lys-C/Trypsin Proteolysis over Trypsin Digestion. J. Proteome Res. 11 (11): 5145-5156 (2012). <https://doi.org/10.1021/pr300273g>)
 - i. Cleaves to the carboxyl side of lysine residues
 - ii. Resistant to chemical denaturation
 - iii. Active at pH 8.0- 9.5
 - iv. Compatible with digestions up to 8 M urea
 - b. Asp-N (Reference: Tarentino, A. L., et al. Molecular Cloning and Sequence Analysis of Flavastacin: An O-Glycosylated Prokaryotic Zinc Metalloendopeptidase. Archives of Biochemistry and Biophysics, 319 (1): 281-285 (1995). <https://doi.org/10.1006/abbi.1995.1293>)
 - i. Cleaves to the amino side of aspartic acid residues
 - ii. Asp-N is a metalloprotease and requires small amounts of zinc to enhance activity
 - iii. Maximal activity at pH 7.0- 8.0
 - c. Glu-C (Reference: Liu, S., et al. Mildly Acidic Conditions Eliminate Deamidation Artifact during Proteolysis: Digestion with Endoprotease Glu-C at pH 4.5. Amino Acids, 48 (4): 1059-1067 (2016). <https://doi.org/10.1007/s00726-015-2166-z>)

- i. Cleaves the carboxyl side of glutamic acid
 - ii. Activity and cleavage specificity is affected by buffer conditions
- d. Chymotrypsin (Reference: Kostka, V. and Carpenter, F. H. Inhibition of Chymotrypsin Activity in Crystalline Trypsin Preparations. JBC, 239 (6): 1799-1803 (1963). [https://doi.org/10.1016/S0021-9258\(18\)91261-5](https://doi.org/10.1016/S0021-9258(18)91261-5))
 - i. Cleaves to the carboxyl side of tyrosine, phenylalanine, tryptophan, and leucine residues
 - ii. Treatment with N-alpha-p-tosyl-L-lysine chloromethyl ketone (TLCK) eliminates residual trypsin activity
- e. PNGase F (Reference: DeRosa, C.M., et al. Simultaneous N-Deglycosylation and Digestion of Complex Samples on S-Traps Enables Efficient Glycosite Hypothesis Generation. ACS Omega 8 (4): 4410-4418 (2023). <https://doi.org/10.1021/acsomega.2c08071>)
 - i. Cleaves to the internal glycoside bond between the asparagine residue and innermost monosaccharide of the N-glycan, N-acetylglucosamine
 - ii. Can hydrolyze high-mannose, hybrid, bi-, tri-, and tetra-antennary oligosaccharides
 - iii. Deamidation modification detected by a mass shift + 0.984 at the deamidation site
 - iv. Active at pH 7.0- 9.0

¹ Reference: Zougman, A., et al. Suspension trapping (STrap) sample preparation method for bottom-up proteomics analysis. Proteomics, 14 (9): 1006-1010 (2014). <https://doi.org/10.1002/pmic.201300553>

² Reference: Chavan, S., et al. Detecting red blood cell protein antigens by tandem mass spectrometry. Transfusion (2025). <https://doi.org/10.1111/trf.18252>

³ Reference: Anmangandla, A., et al. The Acyl-CoA Specificity of Human Lysine Acetyltransferase KAT2A. Biochemistry, 61 (17): 1874-1882 (2022). <https://doi.org/10.1021/acs.biochem.2c00308>

⁴ Reference: Stiving, A.Q., et al. Functionality and translation fidelity characterization of mRNA vaccines using platform based mass spectrometry detection. npj Vaccines 10, 38 (2025). <https://doi.org/10.1038/s41541-025-01082-4>

⁵ Reference: DeRosa, C.M., et al. Simultaneous N-Deglycosylation and Digestion of Complex Samples on S-Traps Enables Efficient Glycosite Hypothesis Generation. ACS Omega 8 (4): 4410-4418 (2023). <https://doi.org/10.1021/acsomega.2c08071>