

## S-Trap™ Micro High Recovery Protocol



- 1) Elute protein from IPs or dissolve protein in 25  $\mu$ L of high recovery urea-SDS lysis/solubilization buffer: 5% SDS, 8 M urea, 100 mM glycine\* pH 7.55.
- 2) Reduce by adding 1  $\mu$ L reductant (120 mM aq. TCEP). Incubate at 37 °C for 15 min.
- 3) Alkylate by adding 1  $\mu$ L alkylator (500 mM MMTS). Incubate at RT for 15 min.
- 4) Add 2.5  $\mu$ L **55% aq. phosphoric acid** to the 25  $\mu$ L sample. This is different than the normal protocol.
- 5) Add 165  $\mu$ L of S-Trap binding buffer (90% MeOH, 100 mM final TEAB, pH 7.55) **into the S-Trap micro column**. It will not flow through.
- 6) **The next two steps must be done as quickly as possible.** Add 2  $\mu$ g of trypsin/lys-C mix into the acidified sample, immediately mix by pipetting up and down, then immediately transfer the mixture into the S-Trap binding buffer held in the micro spin column. Again mix by pipetting up and down.
- 7) Spin in bench-top centrifuge in a standard 1.7 mL sample tube at 4,000 g until all solution has passed through. Remove flow through.
- 8) Wash by adding 150  $\mu$ L S-Trap buffer to the spin column and centrifuging through. Remove flow through. Repeat three times. Protein will not be lost during washes.
- 9) Add 0.5  $\mu$ g of trypsin in 25  $\mu$ L of 50 mM TEAB, pH 8 to the top of the protein trap. The S-Trap is highly hydrophilic and will absorb the solution. However, **ensure there is no bubble atop the protein trap.**
- 10) Cap the spin column **loosely** and incubate in a clean tube for 2 hrs at 47 °C for trypsin. Incubate in a water-saturated atmosphere. **DO NOT SHAKE. The cap MUST NOT form an air-tight seal.**
- 11) Elute peptides with 40  $\mu$ L each of 50 mM TEAB and then 0.2% aqueous formic acid. Add the first TEAB elution to the trypsin solution prior to any centrifugation. Centrifuge elutions through at 4,000 g.
- 12) Elute hydrophobic peptides with 35  $\mu$ L 50% acetonitrile, 0.2% formic acid.
- 13) Lyophilize peptides and resuspend as desired (buffer A or MALDI matrix).

\*Glycine will be removed during washing and helps to limit carbamylation from activated urea.

**PROTIFI**  
Innovative omics solutions